

MANUFACTURE STANDARD, MANUFACTURE RELIANCE

RT-7600S

Auto Hematology Analyzer

NEW



Features


- 20 parameters+3 histograms
- 3-part differentiation of WBC
- 2 counting mode: whole blood and prediluted
- Up to 60 samples per hour
- Color LCD with mouse/keyboard operation
- Built-in thermal printer, external printer optional
- External AC adapter avoids electronic noises, enhancing accuracy
- Up to 35000 sample results (including histograms) can be stored
- Multi-language software available on request
- Original factory's blood control and calibrator



website: www.medi-world.net

RT-7600S

Auto Hematology Analyzer



Technical Specifications

Parameters	WBC, LYM#, MID#, GRA#, LYM%, MID%, GRA%, RBC, HGB, MCHC, MCH, MCV, RDW-CV, RDW-SD, HSB, LAN, keyboard and mouse interface
Principles	Electrical resistance for counting WBC, RBC and PLT
Sample Volume	Prediluted 20µL, Whole Blood 9.8µL
Throughput	Up to 60 samples per hour
Display	8" color LCD
Alarms	Error messages
Input/Output	RS-232, USB, LAN, keyboard and mouse interface
Printout	Built-in thermal printer, external printer optional
Power Requirement	a.c.110V-220V, 50Hz/60Hz
Operating Environment	Temperature: 15°C~35°C; Humidity: 10%~90%
Net Weight	18kg
Dimension L×W×H (mm)	436×363×367



Rayto Life and Analytical Sciences Co.,Ltd.

2nd Floor, 5th Building, Software Park, NO.2, Gaoxinzhong
3rd Road, High-tech Industrial Park, Nanshan,
Shenzhen 518057, P.R.China
Tel:+86(755)86868996 Fax: +86(755)86168796
Email:info@rayto.com

