



www.biocare.com.cn



ECG-1210

Digital Electrocardiograph



ECG-1210

12-channel digital electrocardiograph


● Data management solution

- Internal memory can store 3000 files
- Support SD card and USB flash disk
- Upload data to ECG management software via WIFI (optional)





- Upload data to ECG management software by LAN or WIFI (optional)

● HL7 for seamless data transfer



● Streamlined workflow from start to finish

- Icon buttons on the screen for intuitive operation
- Various report formats for printing



- Barcode scanner to quickly input patient ID

- On-screen guide to check for loose leads

- Preview the report before printing

Shenzhen Biocare Bio-Medical Equipment Co., Ltd.

Add: #16-1, Jinhui Road, Jinsha Community, Kengzi Sub-District, Pingshan New District, 518122 Shenzhen, P.R. CHINA

Tel: +86-755-2796 0888/3300 5899

Fax: +86-755-2796 0643

E-mail: sales@biocare.com.cn

Website: www.biocare.com.cn

© 2019 BIOCARE all rights reserved. Specifications subject to changes without prior notice.

ENG-CATA-V2.6-20190509